



margin medical instruments

*ENT treatment system*

## ***MT 1500***

The MediTop 1500 ENT Treatment system has been specially designed for Ear Nose and Throat departments out patient facilities. It's 'Space Saving Concept, allows rooms to be utilised far more efficiently. The 1500 System can be custom made to your own exact specifications and a room layout design will be done to ensure the 1500 System provides an efficient working environment for Consultants and Patients. The 1500 System is also available as a mobile system.








### Features:

- Two sliding instrument trays.
- One twin lower used instruments and waste disposal container featuring an open chute for easy access.
- Large storage drawer.
- Mains supply 230 Volts - 50 Hz
- Two 230 Volts outlet sockets.
- In standard format the unit will be either grey or white with lacquered front panels in a range of five different RAL colours.



Main control display featuring nine function switches and four instrumentholders (three of them switchable)

	RAL 5015
	RAL 6027
	RAL 5021
	RAL 9003
	RAL 7035

Colour can slightly differ because of printing process. On request other colours can be delivered against an extra charge.



### Options:

- Desktop in beech wood laminate;
- One large instrument tablet instead of the standard subdivision;

### Built in options:

- Suction system;
- Compressed air system;
- Mirrorheater;
- Headlight supply;
- Switchable hook;
- Fibre optic system;
- Warm water syringe;
- Pre-Mirrorheater;
- Cautery system;
- Scope holders.
- Holder for built in microscope.

### Dimensions:



margen medical instruments

Crowthorne Enterprise Centre, Old Wokingham Road,  
Crowthorne, Berkshire RG45 6AW  
Tel: 01344 751678 Fax: 01344 751601  
Email: [info@margenmedical.co.uk](mailto:info@margenmedical.co.uk)  
[www.margenmedical.co.uk](http://www.margenmedical.co.uk)