



margin medical instruments


ENT treatment system

MT 750

The space saving design concept of this mobile ENT treatment system has been achieved through years of research and development.

It improves the working environment for the ENT consultant and nursing staff. Also increased efficiency for patients in busy outpatients clinics.

In short this system will improve the work environment for staff and less waiting time for patients.



Features:

- One sliding instrument tray;
- One fixed instrument tray;
- One lower used instruments container and waste disposal bucket featuring an open chute for easy access;
- Large storage drawer;
- Two function outlet sockets;
- Standard unit colour white or grey
- Optional coloured front panels lacquer finished




Built in options:

- Suction system;
- Compressed air system;
- Mirrorheater;
- Headlight supply;
- Switchable hook;
- Fibre optic system;
- Warm water syringe;
- Pre-Mirrorheater;
- Cautery system;
- Scope holders.



Main control panel featuring four instrumentholders. (three of them switchable)

Dimensions:



margen medical instruments

Crowthorne Enterprise Centre, Old Wokingham Road,
Crowthorne, Berkshire RG45 6AW
Tel: 01344 751678 Fax: 01344 751601
Email: info@margenmedical.co.uk
www.margenmedical.co.uk